

**AUTOMATIC BIO MEDICAL
LIQUID WASTE TREATMENT SYSTEM**



Regular Maintenance Visit



End User Training



Installation Done by
Trained Professionals



Compact Equipment



Touch screen LCD display



- Fully Automatic Liquid Waste Treatment System with manual operating mode.
- Continuous treatment cycle to ensure uninterrupted working with an option to treat any level of liquid waste manually.
- PLC based touch screen LCD display to operate equipment.
- Functions like starting the cycle, dispensing sodium hypochlorite solution, drainage etc. are automatically performed by the equipment.
- Data logging using USB 2.0.
- The equipment gives print ticket with cycle data for record keeping.
- User login is protected by a pass code.
- pH Level monitoring.
- Audio visual alarms for low / full levels of hypochlorite solution, full level for treatment and storage tanks.
- Tested for efficacy by Shriram Institute for Industrial Research, New Delhi with 99.99% efficacy.

Model	Capacity
Hypogermished 150	150 Litres



The Equipment treats and disposes of hospital liquid waste as per Biomedical Waste Management Rules 2016 and CPCB Guidelines.